

Day 2 - LoDoPaB Dataset Exploration

In this notebook, we will:

1. Load the **LoDoPaB** dataset using `dival`.
2. Inspect the dataset statistics (train/validation/test sizes).
3. Sample and visualize data (ground truth vs. sinogram).
4. Apply HU clipping and normalization as preprocessing.

1. Import packages & load dataset

We use `LoDoPaBDataset` from `dival.datasets`. Since ASTRA CUDA backend is not available in this environment, we explicitly set the implementation to `skimage` (CPU).

```
In [3]: from dival.datasets import LoDoPaBDataset
import matplotlib.pyplot as plt
import numpy as np
import random
```

```
ds = LoDoPaBDataset(impl='skimage')
```

The LoDoPaB-CT dataset could not be found under the configured path 'D:\cosc_4372\projects\lowdose_ct_project\ct-reconstruction-pipeline\data\raw\lodopab\y'.

Do you want to download it now? (y: download, n: input other path)

Path to LoDoPaB dataset:

updated configuration in 'C:\Users\micke\.dival\config.json':

'lodopab_dataset/data_path' = D:\cosc_4372\projects\lowdose_ct_project\ct-reconstruction-pipeline\data\raw\lodopab

2. Dataset statistics

Check the number of samples available in each split (train/validation/test). This ensures the dataset has been loaded correctly from the configured path.

```
In [4]: print("Train:", ds.get_len("train"))
print("Validation:", ds.get_len("validation"))
print("Test:", ds.get_len("test"))
```

Train: 35820

Validation: 3522

Test: 3553

3. Inspect a sample

We retrieve a sample from the training set and explicitly convert it to `numpy.ndarray`. This allows us to examine array shapes and intensity ranges for both ground truth (CT slice) and sinogram (projection data).

```
In [5]: observation, ground_truth = ds.get_sample(0, part='train')

to_np = lambda x: x.asarray() if hasattr(x, "asarray") else np.asarray(x)
sino_np = to_np(observation)
gt_np = to_np(ground_truth)

print(type(gt_np), gt_np.shape, gt_np.min(), gt_np.max())
print(type(sino_np), sino_np.shape, sino_np.min(), sino_np.max())

<class 'numpy.ndarray'> (362, 362) 0.0 0.489164
<class 'numpy.ndarray'> (1000, 513) -0.000798698 0.0802129
```

4. Visualization

Display the **Ground Truth CT slice** and the corresponding **Low-dose Sinogram** side by side. This step provides an intuitive understanding of how the data looks before preprocessing or reconstruction.


```
In [7]: plt.figure(figsize=(8,3))

plt.subplot(1, 2, 1)
plt.imshow(gt_np, cmap='gray')
plt.title('Ground Truth')
_ = plt.axis('off')


plt.subplot(1, 2, 2)
plt.imshow(sino_np, cmap='gray', aspect='auto')
plt.title('Sinogram')
_ = plt.axis('off')

plt.tight_layout()
plt.show()
```

Ground Truth



Sinogram



```
In [9]: for i in range(5):
        idx = random.randint(0, ds.get_len("train")-1)
        obs, gt = ds.get_sample(idx, part='train')
        obs_np = obs.asarray() if hasattr(obs, 'asarray') else np.asarray(obs)
        gt_np = gt.asarray() if hasattr(gt, 'asarray') else np.asarray(gt)


        plt.figure(figsize=(8, 3))

        #Plot ground truth
        plt.subplot(1, 2, 1)
        plt.imshow(gt_np, cmap='gray')
        plt.title(f'Ground Truth {idx}')
        _ = plt.axis('off')


        #Plot Sinogram
        plt.subplot(1, 2, 2)
        plt.imshow(obs_np, cmap='gray', aspect='auto')
        plt.title(f'Sinogram {idx}')
        _ = plt.axis('off')

        plt.tight_layout()
        plt.show()
```


Ground Truth 16558




Sinogram 16558




Ground Truth 29243




Sinogram 29243




Ground Truth 2011




Sinogram 2011




Ground Truth 29754



Sinogram 29754



Ground Truth 690



Sinogram 690

